

Magnet Sales Information


500 MHz / 54 mm

ASCEND™



Magnet System 500/54 ASCEND™
Part Number Z115311

Geometrical Dimensions



⊙ = magnetic center
1 mT = 10 Gauss

Z1044563



Geometrical Dimensions

		Description
C =	2520 mm	Operational ceiling height
D-RT =	54 mm	Diameter room temperature bore tube
D1 =	745 mm	Diameter cryostat upper part
D2 =	795 mm	Diameter cryostat bottom plate
H1 =	1005 mm	Height of cryostat from bottom flange – upper flange
H2 =	1391 mm	Height of cryostat from bottom flange to Helium tower Minimum height for transportation
H3 =	1564 mm	Height of cryostat from bottom flange to Helium manifold
S =	570 mm	Height between floor and magnet bottom flange


System Data

Minimum operational ceiling height (Helium transfer line 29085)	2520 mm
Minimum ceiling height with standard Helium transfer line 53962	2805 mm
Required space (footprint)	~ 1.4 m ²
System weight (empty, without magnet stand)	440 kg
Magnet stand easy type / Optional damping cylinders	32 kg / 120 kg
System weight (filled completely, with magnet stand)	650 kg / 740 kg

NMR Magnet Specifications

Type	BZH 500'70 ASCEND™
NMR-frequency (¹ H)	500 MHz
Operating field	11.74 Tesla
Field stability (guaranteed value in persistent mode)	< 10 ppb/hr (< 5.0 Hz/hr)
Axial range with homogeneity better than 10ppm	~ 55 mm
Radial stray field (horizontal distance of the 0.5mT (5G) line from the magnetic centre)	< 0.60 m
Axial stray field (vertical distance of the 0.5mT (5G) line from the magnetic centre)	< 1.20 m
Cryo shims	X, Y, Z, Z ² , Z ³ , XZ, YZ, X-Y, X ² -Y ²
Electromagnetic Disturbance Suppression EDS** typical	> 99 %

Fringe Field Plot



Fringe field contour	Radial [R]	Axial [Z]
200 mT (Directive 2004/40/EC)	Inside cryostat	Inside cryostat
5 mT	0.47 m	0.80 m
3 mT	0.50 m	0.87 m
1 mT	0.56 m	1.06 m
0.5 mT (5 Gauss)	0.60 m	1.20 m
0.2 mT	0.74 m	1.48 m
0.1 mT	0.96 m	1.75 m
0.05 mT (~Earth magnetic field)	1.26 m	2.09 m



Cryostat Specifications


Type	D 315/54 ASCEND™
Room temperature bore	54 mm
Approx. Helium evaporation rate under stabilized conditions (T=20°C, p=1030 mbar)	<13 ml liquid Helium/hour
Liquid Helium refill volume/total volume	~ 57 /85 litres
Helium hold time	> 180 days
Approx. Nitrogen evaporation rate under stabilized conditions (T=20°C, p=1030 mbar)	<230 ml liquid Nitrogen/hour
Liquid Nitrogen refill volume/total volume	~ 85/106 litres
Nitrogen hold time	> 15 days

Accessories

Magnet stand H80 570-720 EMI Easy Easy installation of magnet system with Elastomeric Isolators Frequencies damped > 14 Hz / resonance frequency = 9.5 Hz	Standard	Z123639
Magnet stand F80-570 ADI (height 570 mm) Air Spring Damped Isolators with vertical damping Frequencies damped > 3.8 Hz / resonance frequency = 2.6 Hz	Optional AH0063	Z112742
Magnet stand F80-700 ADI (height 700 mm) Air Spring Damped Isolators with vertical damping Frequencies damped > 3.8 Hz / resonance frequency = 2.6 Hz	Optional if S2 Probes are used	Z112741
Magnet stand F80-570 API (height 570 mm) Air Piston Damped Isolators with vertical and horizontal damping Frequencies damped > 3.8 Hz / resonance frequency = 2.6 Hz	Optional	Z117691
Magnet stand F80-700 API (height 700 mm) Air Piston Damped Isolators with vertical and horizontal damping Frequencies damped > 3.8 Hz / resonance frequency = 2.6 Hz	Optional if S2 Probes are used	Z117692
API cylinders upgrade when magnet stand F80 is existing Air Piston Damped Isolators with vertical and horizontal damping Frequencies damped > 3.8 Hz / resonance frequency = 2.6 Hz	Optional AH0065	Z114651
Electronic atmospheric pressure device with adjustable set point*		Z102597
Nitrogen level sensor for BSMSII (SCB3) bendable (not good for BSNL)	Included	Z122394

Equipment for Cryogen Transfer

Helium transfer line* D3xx (1455/2060/655)	AH0070	53962
Helium transfer line* with bendable extensions (1455/2060/380) for minimum operational ceiling height (2850 mm)		29085



Geometrical Dimensions with F80 Magnet Stand

Width of magnet stand	W1	1236 mm
Depth of magnet stand	W2	1116 mm
Diameter of magnet stand = 2 x radius	D3	1375 mm

Transport

Overall system dimensions for transportation		
Magnet box	L x D x H	114 x 93 x 198 cm ³
Magnet stand box (Optional LC Stand)	L x D x H	97 x 76 x 120 cm ³
Magnet stand box (Standard stand)	L x D x H	125 x 77 x 60 cm ³
Minimum system dimensions of magnet, unpacked (without manifold)		Ø 79.5 cm, H 139 cm
Magnet System weight for transportation		~ 670 kg
Magnet stand box weight for transportation (standard stand / optional)		~ 80 / 160 kg

Installation

Liquid Nitrogen needed for cooling down	400 litres
Liquid Helium needed for cooling down	200 litres
Liquid Helium needed for energizing, cryo shimming and quench reserve	300 litres
Nitrogen gas for flushing	1 cylinder 50 l/200 bar
Helium gas for flushing	1 cylinder 50 l/200 bar

* A detailed description of the marked objects can be found in chapter „Accessories“.

Electromagnetic Disturbance Suppression **EDS™:

Ascend™ magnets efficiently suppress external electromagnetic field disturbances using a proprietary Bruker technology **EDS™**. Sources of such disturbances are corridor traffic, elevators, power lines, outside vehicular traffic and railway lines.

Definition:

The **EDS™ factor** for spatially homogeneous disturbances is defined as the fraction of the external disturbance suppressed by the magnet in the magnetic centre at a given disturbance frequency. Thereby, no digital lock system or other field compensation device is used.

Detailed specification for magnet system 500/54 Ascend™:

Disturbance frequency	EDS™ factor
< 0.01 Hz	> 99 %
0.01 - 1 Hz	> 97 %
1 - 5 Hz	> 98 %
> 5 Hz	> 99 %
16.667 Hz (railways)	> 99.9 %
50 Hz (power lines)	> 99.9 %
60 Hz (power lines)	> 99.9 %





- **Bruker BioSpin**
your solution partner

Bruker BioSpin provides a world class, market-leading range of analysis solutions for your life and materials science needs

- **Bruker BioSpin Group**


sales@bruker-biospin.com
www.bruker-biospin.com

LISTA CZĘŚCI		
ELEMENT	NUMER KATALOGOWY	OPIS
1	65018687	BIS RapidStrut 41x82x2,5 BUP
2	65018647	BIS RapidStrut 41x41x2.5 BUP
3	665885400	BIS RapidStrut® Uchwyty ściennie G2 BUP
4	665885408	BIS RapidStrut łącznik siodłowy 41x82
5	6621083	BIS Wspornik do konsol ściennych regulowany
6	61481030	BIS Śruba M10x30
7	65381010	BIS podkładka M10 BUP
8	65188010	BIS Nakrętka ślizgowa nie zmont. M10



UWAGI:

1. Kotwić do ściany za pomocą kotwy chemicznej WVSF200 (nr. kat. 6099125E) i trzpieni M10 (nr. kat. 60991013).

Opracował: mgr inż. Grzegorz Sobczyk	Sprawił: Przed montażem należy sprawdzić rysunek i zweryfikować z ofertą.	Date 04.06.2020	Skala: -	Format: A3
			Tytuł: Stelaż centrali mocowany do ściany	
			Nr rys: 200604-1	
Objekt: -			Arkusze 1 / 1	

Niniejszy rysunek stanowi informację o produktach firmy Walraven Sp. z o.o. i nie stanowi projektu w rozumieniu właściwych przepisów. Opracowanie stanowi własność Walraven Sp. z o.o. Montaż nie zgodny z wytycznymi na rysunku, nie gwarantuje prawidłowego funkcjonowania. Walraven nie ponosi odpowiedzialności za jakiegokolwiek szkody wynikające z doradztwa lub informacji udzielonych jej klientom, o ile wyraźnie pisemnie nie uzgodniono inaczej.